MÉRÉSTECHNIKA

MEASURING TECHNICS

To Metrology and Beyond

New technical advances take industrial surface inspection to the next level

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Keywords: laser scanning confocal microscope, metrology, resolution, magnification

Introduction

In the modern industrial era, where perfection is expected by the customer, manufacturers are facing increasing pressure to produce components of the highest quality. This coupled with the push for miniaturisation and complexity adds pressure on the R&D and quality checking stages. New product output can therefore be slowed down as researchers check every novel component carefully before going to production. Traditionally, these microscopic features are viewed on instruments which leave a sample unusable after testing; the scanning electron microscope (SEM) for instance. Metrology requires the use of more than just an SEM though, with stereo microscopes and contact profilometers also featuring heavily. To eliminate the need for multiple instruments and destructive processes, a new metrology instrument has been developed by Olympus, which combines a conventional light microscope and a laser scanning confocal microscope.

Analyse, don't annihilate

Metrological analyses are best carried out using the finest possible resolution at a range of magnifications. At present, for many high magnification observations samples are viewed using an SEM. This is a time consuming and destructive process, since the sample must be broken up, coated with a conductive material and can only happen under a vacuum. SEM images also provide very limited opportunities for analysis, acting more as an observational tool than a measuring device. Before reaching the SEM stage, the sample may be viewed under a conventional stereo microscope to look for more obvious surface details. This microscope also acts as a useful tool during the preparation of a sample for SEM or contact profilometry.

OLYMPUS LIFE AND MATERIAL SCIENCE EUROPA GMBH
Microscopy Esther Ahrent Section Manager Marketing
Communication Tel: +49 40 2 37 73 – 5426 Fax: +49 40 2 37
73 – 647 E-mail: microscopy@olympus-europa.com

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The LEXT is designed to essentially replace these three different tools with one single solution. This has required the application of a number of technical advances which allow the system to move easily from lowest magnification – 120x, to the highest 14,400x, whilst providing peerless resolution and analysis capabilities.

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Figure 1: Olympus LEXT – High resolution 3D surface profile observation and measurement with UV-laser confocal microscopy

Larger & Clearer

At 120x the LEXT system overlaps with the magnification provided by stereo microscopes making it easy to carry out initial observation of conspicuous surface details. Areas of interest can then be viewed at much higher magnifications (up to 14,400x) which overlap with SEM, by moving through a series of automated parfocal objectives. These are designed so that when changing magnification, no re-focusing is required. This means that keeping track of a specific feature is much easier than normal and also that the user can get as close as they want to the sample without breaking it apart or causing any other damage.

Standard laser confocal systems can achieve resolutions down to 200nm, but the technological advances in the LEXT allow it to reach as low as 120nm in x-/y-direction. Moreover, resolution on the Z-axis is 10nm, due to the use of a new mathematical algorithm and the finely controlled movement of the nosepiece rather than the object stage. In addition, the LEXT used a 408nm laser diode with newly designed and calibrated optical components, which have excellent transmission across the entire visible spectrum.

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
Visualisation

By covering such a wide magnification range with parfocal objectives and combining a conventional and laser confocal microscope, viewing a sample on this system is very easy. The stage can handle items up to 10kg in weight and 100mm in height, and therefore a complete brake disc or turbine blade can be analysed with the same ease as a micro-electromechanical system (MEMS) device. Using conventional microscopy specific features of a surface can be quickly identified at low magnification. The areas of interest can then be magnified further and analysed using a plethora of powerful, confocal microscopy-based observation and measurement tools.


With the majority of surfaces being analysed for contours, special techniques need to be employed at high magnifications since there is very little depth of focus. The LEXT uses confocal optics, each optical slice consists only of in-focus pixels and therefore, by simply combining a stack into one image, a perfect 2D representation of the entire feature is formed. For many of the analyses carried out using LEXT though, it is important to see the object in 3D detail and therefore these same optical slices can easily be arranged into a highly accurate 3D image. Most confocal microscopes carry out these processes post-capture and therefore there is a significant delay between taking the images and seeing the result. Once the top and bottom of the stack has been set the LEXT does not have this delay since it uses a specialised MEMS-based scanning unit and advanced software providing real-time observation.

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Figure 3: Ultra precise drill, depth and angle measurement

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
Figure 2: Combined 408 nm laser confocal and microscope optical system

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Figure 4: IC Stud Bump:

a) real colour 3D image; b) 3D wireframe image with colour coded height information; C) 3D greyscale image

Analysis and reporting


The images produced can be shown in a number of different ways. Not only the standard 2D and 3D images mentioned above, but also as wire mesh and colour contour images. The user can then select areas or cross sections to analyse for height, length, angle and volume data as well as surface roughness. Another impressive feature is the ability to combine the conventional microscopy channel output with that of the laser confocal; this means that true colour can be used to show clearer contrast. This is especially useful if a sample consists of a number of areas at equal heights, but made from different materials. This type of analysis cannot be completed on SEM since the surfaces are all coated with the same compound and therefore would appear to have the same height.

With repeatability and traceability key to the auditing of most manufacturing and R&D procedures, it is essential that any measuring instrument can store the correct information about each and every procedure. This must also be available in printed reports. All analyses completed by the LEXT system conform to international standards such as the United Kingdom Accreditation Service (UKAS), the Japan Quality Assurance Organisation (JQA) and the Physikalisch-Technische Bundesanstalt (PTB). Bespoke standards can also be programmed into the system without limiting its operability. Comprehensive reporting via a built-in facility provides users with the flexibility to include all of the data required in a consistent format.

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Conclusions

Whether analysing the roughness of a vehicle brake disc or the angular intricacies of a MEMS micro-weld, the observations required should not cause any damage to the sample – especially during R&D stages where there is often only one unit. Many current procedures though, do exactly this which can escalate the development and production costs. The application of unique optics and specialised mathematical algorithms is set to change all this in the new LEXT from Olympus. Full metrological analysis may currently require the use of a stereo microscope, SEM and contact profilometer. The LEXT can essentially do the work of all three instruments since it covers a magnification range which overlaps with both stereo and scanning electron microscopes. It is also capable of very comprehensive analysis conforming to all international standards and replaces the need for a contact profilometer.

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